TECHNOJET

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A DURABLE JET-NEBULIZER WITH PARTICLE SELECTOR AND VALVE SYSTEM

TECHNOJET IS A NEWLY-DESIGNED JET NEBULIZER CAPABLE OF DELIVERING ANY MEDICATION IN A FAST AND EFFICIENT WAY. THE SELECTOR CONSENTS TO CHANGE PARTICLE DIMENSION REACHING DIFFERENT AREAS OF THE RESPIRATORY TRACT AS FOLLOWS:

POS. 1: PARTICLES WITH DIAMETER BETWEEN 1 AND 5 MICRON REACH THE ALVEOLI.

POS. 2: PARTICLES WITH DIAMETER BETWEEN 6 AND 9 MICRON REACH LOWER AIRWAYS (TRACHEOBRONCHIAL TRACT).

POS. 3: PARTICLES WITH DIAMETER OVER 10 MICRON REACH UPPER AIRWAYS (NOSE, LARYNX AND TRACHEA).


THE SELECTOR PERMITS TO CHOOSE PARTICLE DIMENSION WITH AN EASY MOVE, AVOIDING INTERCHANGING PARTS. THANKS TO THE REDUCTION OF WASTED MEDICATION, TECHNOJET IS ALSO RECOMMENDED FOR SEVERE RESPIRATORY DISORDERS SUCH AS COPD OR ASTHMA, WHEN EXPENSIVE MEDICATIONS ARE PRESCRIBED OR FOR INTENSIVE USE IN GENERAL.

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
TO APPRECIATE THOROUGHLY THE PERFORMANCES AND BENEFITS OF 3A NEBULIZERS IT IS ESSENTIAL TO KNOW THE FOLLOWING CONCEPTS OF AEROSOL THERAPY:

CONCEPT	MEANING	EXPRESSED
Neb rate	Amount of aerosol delivered for minute	ml/min
Respirable fraction	Amount of particles below 5 micron and therefore capable of reaching deep airways	% particles
MMAD o MMD	Mass Median (Aerodynamic) Diameter. It indicates that 50% of the total volume nebulized is in particles with diameter lower than the MM(A)D	micron μ m
Residual volume	Amount of medication left in the nebulizer. It cannot be nebulized.	ml

- POS. 1
 NEB RATE: 0.20 ML/MIN
 PARTICLE SIZE (MMAD): $1 \pm 5 \mu\text{m}$
- POS. 2
 NEB RATE: 0.35 ML/MIN
 PARTICLE SIZE (MMAD): $6 \pm 9 \mu\text{m}$
- POS. 3
 NEB RATE: 0.50 ML/MIN
 PARTICLE SIZE (MMAD): $10 \pm 14 \mu\text{m}$

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NEBULIZER TECHNICAL DATA

CAPACITY	12 ml
RESPIRABLE VOLUME $< 5 \mu\text{m}$	80%
RESIDUAL VOLUME	0.7 ml

CE 0051
(dir. 93/42/EEC)

TECHNOJET